## ISO Cube<sup>TM</sup> Daily QA Package

Model 023-05



### AFFORDABLE "TURNKEY" SOLUTION FOR DAILY MACHINE QA

Target positioning through imaging guidance is critical for the accurate delivery of radiation treatment. Verifying that all of the imaging, localization and targeting systems are aligned with the true radiation isocenter is crucial. The CIRS ISO Cube™ Daily QA Package provides a cost-effective, fast and accurate means of testing radiation isocenter coincidence with the isocenters of the image guidance systems.

The package includes ISO Cube<sup>™</sup> Daily QA Phantom, ISO Base<sup>™</sup> Platform and ISO Analyze<sup>™</sup> Image Analysis Software.

ISO Analyze<sup>™</sup> integrates with the ISO Cube<sup>™</sup> and ISO Base<sup>™</sup>, enabling user-friendly quality control of the isocenter of a LINAC by analyzing DICOM images acquired with the EPID detector. Controls are run automatically, analyzing images of the phantom and quantifying a large number of evaluation parameters. The software allows users to easily generate, save and print a report for each preceding control.

The ISO Base<sup>™</sup> is used to position and level the ISO Cube<sup>™</sup> on the treatment couch. It contains integrated pixel calibration targets for use with ISO Analyze<sup>™</sup>.

2428 Almeda Avenue Suite 316 • Norfolk, Virginia 23513 • USA Tel: 800.617.1177 • 757.855.2765 • Fax: 757.857.0523

WWW.CIRSINC.COM

The phantom, base and software were designed specifically for daily system checks. LINAC laser and light fields can be "tuned" to true radiation isocenter using the engraved markings on the exterior of the phantom. The light field and radiation field alignment can be checked using the phantom's integral radiographic markers. More importantly the isocenters of both the OBI and the EPID can be checked for true special alignment and coincidence with that of the treatment beam.

The ISO Cube contains a center target and an offset target. Both targets measure 6.35 mm in diameter and are made of ceramic. The center target imaged in concert with the external concentric engraved circles provides greater accuracy localizing the center of the phantom with respect to the center of the radiation field. The offset target is used to insure the table offset coordinates generated by kV/MV imaging are accurate.



# ) **ISO** ) ANALYZE

POWERED BY AQUILAB

# Integrated Data Analysis in 5 Easy Steps

#### **ACQUIRE**

Use pre-established ISO Cube's treatment plan to acquire all necessary images for analysis of parameters.

ISO Analyze<sup>™</sup> Software provides suggested image sequences for each parameter









#### DEFINE

One time setup for each LINAC requires:

- LINAC Identity
- Acquisition Conditions
- Acquisition Sequences

![](_page_1_Picture_14.jpeg)

Name*	Norfolk General
Manufacturer*	Varian
Model*	Trilogy
Serial number*	XXX

![](_page_1_Picture_16.jpeg)

#### SELECT

Select pre-acquired DICOM images to use in each parameter calculation

Select images are highlighted for each parameter

![](_page_1_Picture_20.jpeg)

### WHY USE ISO CUBE<sup>™</sup>?

- 1. The stereo-triangulation approach employed with ISO Cube is similar to that of highly accurate GPS Systems. The calculation accuracy of the isocenters is based on projections of targets distributed in a 3D space contained within the 3D space found at the convergence of the treatment beams.
- 2. ISO Cube allows for comprehensive and thorough testing of LINAC specific centers of rotation.
- 3. ISO Cube allows for quick, daily QA of OBI and CBCT isocenters.
- ISO Cube shows the impact of posi-4. tioning a 3D object on the couch using LINAC specific positioning lights and lasers.
- 5. ISO Cube offers an opportunity to assess the orthogonality of LINAC specific positioning lights and lasers.
- 6. ISO Cube allows calculation of the planar deviation between the center of the radiation field and the projection of the LINAC mechanical ISO center on the image plane for each gantry and collimator position. This eases troubleshooting of clear outliers.
- 7. ISO Cube's internal features permit assessment of misalignment between light field and radiation field.

#### CALCULATE

ISO Analyze<sup>™</sup> quickly calculates parameters and displays detailed results

![](_page_1_Picture_31.jpeg)

![](_page_1_Figure_32.jpeg)

![](_page_1_Picture_33.jpeg)

![](_page_1_Picture_34.jpeg)

![](_page_1_Picture_36.jpeg)

#### REPORT

Detailed results can be manipulated for indepth analysis

Summary report is displayed and PDF document may be printed and saved

	<b>RS</b> ⊕	ISO ISO Date of ca	Analyze Repor	t
	POWE	RED BY AQUILAB		
Tested device	•			
Name	Norfolk General			
Manufacturer	Varian			
Model	Trilogy			
Serial number	3000			
Phantom info	rmation			
Phantom serial	XXX			
Central marker				
type	Ceramic marker			
Pixel size of t	he MV detector	Computed Value	Theoretical Value	
X pixel size (mm)		0.5292	0.5227	
Y pixel size (mm)		0.5285	0.5227	
Y pixel size (mm) Pixel size of ti Pixel size of ti	of the kV detector	0.5285	0.5227	
Y pixel size (mm) Pixel size of the pixel size size of the pixel size size size size size size size size	of the kV detector	0.5385	0.5227 Theoretical Value	
Y pixel size (mm) Pixel size of t Pixel size of t Pixel size in u (mm)	f the kV detecto	0.5285 Computed Value 0.2619	0.5227 Theoretical Value 0.2586	
Y pixel size (mm) <b>Pixel size of</b> Pixel size of the Pixel size in u (mm) Pixel size in v(mm)	of the kV detector	0.5385	0.5227 Theoretical Value 0.2586 0.2586	
Y pixel size (mm) Pixel size of ti Pixel size is u (mm) Pixel size is u (mm) Collimator Collimator rot	of the kV detector	0.5385  Competed Value  0.2699  0.362  2  2  T  Tel Fiducial/Marker	0.5227 Theoretical Value 0.2566 0.2566	
Y pixel size (mm) Pixel size of the pixel size of the pixel size of the pixel size is u (mm) Pixel size is u(mm) Collimator rot Contract Fukusia/Parker	f the kV detector	0.5385  Competed Value 0.5459 0.5459 2  F  F  Competed Value 0.5459 0.5452 2  C  C  C  C  C  C  C  C  C  C  C  C	0.5227 Theoretical Value 0.2565 0.2565	
Y pixel size (mm) Pixel size of ti Pixel size of ti Pixel size is u (mm) Pixel size is u (mm) Pixel size is u(mm) Collimator rot Collinator Visionia/Marker X (mm)	f the kV detector	0.5385 r Competed Value 0.2489 0.2489 0.2622 2 r al Fiducial/Marker 0.0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	0.5227 Theoretical Value 0.2566 0.2566	
Y pixel size (mm)  Pixel size of ti  Pixel size of ti  Pixel size is u (mm)  Pixel size is u (mm)  Pixel size is u (mm)  Collimator rot  Collimator rot  Contral FiducialNafare  Y (mm)  Contral FiducialNafare	f the kV detector	0.5385  Competed Value 0.2699 0.582  27  ral Fiducial/Marker 0.4  0.4  0.4  0.4  0.4  0.4  0.4  0.4	0.5227 Theoretical Value 0.2586 0.2586 771 771 771	
Y pixel size (mm) Pixel size of ti Pixel size of ti Pixel size is (mm) Pixel size is (mm) Collimator rot Contact Faksik/Narker Conta	of the kV detector	0.5385 7 Competed Value 0.5489 0.5489 0.5422 2 F F Fiducial/Marker 0.0 0 0 0 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	0.5227 Theoretical Value 0.2566 0.2566 7711 227 246 703	
Y plead size (wm)  Plycel size of U  Proved size is (wm)  Proved size is (wm)  Collimator rot Collimator Col	f the kV detector	0.5385 r Competed Value 0.2489 0.2627 27 rral Fiducial/Marker 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0.0 0.	0.5227 Theoretical Value 0.2566 0.2566 7711 227 766 707 797 709 797 60 701 797 60 701 797 60 701 797 60 701 797 60 701 797 60 701 797 60 701 797 60 701 797 60 701 797 60 701 797 60 701 707 70 70 70 70 70 70 70 70 70 70 70 70	
V pland size (mm)  Plykel size of fil  Plykel size of fil  Pland size is (mm)  Pland size is (mm)  Collimator Collimator fol  Collimator fol  Collimator fol  Collimator fol  Collimator fol  Beam Pland Size Size  Common Consult Plands  Demons Consult Pl	f the kV detector	0.5385  Compared Value 0.269 0.262 2  F  ral Fiducial/Marker 0.0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	0.5227 Theoretical Value 0.2566 0.2566 2.256 237 246 257 257 257 257 257 257 257 257 257 257	

#### ISO CUBE™ DAILY QA PACKAGE

#### Model 023-05

![](_page_2_Figure_2.jpeg)

Figure 1. 2D/2D match of kV and DRR.

![](_page_2_Picture_4.jpeg)

Figure 2. Concentric circles verify accurate alignment of ISO Cube and establish true position of the kV radiation isocenter.

#### **SPECIFICATIONS**

OVERALL DIMENSIONS:	12 cm x 12 cm x 12 cm (4.75" x 4.75" x 4.75")
WEIGHT:	1.7 kg (3.7 lb)
MATERIALS:	Plastic Water®

#### **MODEL 023-05 INCLUDES**

QTY	DESCRIPTION
1	ISO Cube™ Daily QA Phantom
1	ISO Base™ Alignment Platform
1	ISO Analyze™ Image Analysis Software
-	User Guide
-	48 Month Warranty

#### **OPTIONAL ACCESSORY**

PART NO.	DESCRIPTION
023-01	Iso Opt Stereotactic Target Frame Adapter

![](_page_2_Figure_12.jpeg)

![](_page_2_Figure_13.jpeg)

![](_page_2_Figure_14.jpeg)

![](_page_2_Figure_15.jpeg)

![](_page_2_Figure_16.jpeg)

©2013 Computerized Imaging Reference Systems, Inc. All rights reserved. All brand names, product names or trademarks belong to their respective holders. Specifications subject to change without notice. Publication: 023-05 DS 091913

![](_page_2_Picture_18.jpeg)

Computerized Imaging Reference Systems, Inc. has been certified by UL DQS Inc. to (ISO) 9001:2008. Certificate Registration No.10000905-QM08.