


Which C-Arm Table is right for you?

Pain Management C-Arm Table - 870

Ideal for pain care applications



SPECIFICATIONS:

Dimensions: 93" l x 26" w (233 x 66 cm) with OR accessory rails
 Tabletop: 90" l x 24" w (228.6 x 61 cm)
 Radiolucent Area:
 60" l x 24" w (152.4 x 61 cm) unobstructed head to accessory rail
 4" l x 17.5" w (10.1 x 44 cm) additional, between accessory rails
 Radiolucent Area, Total Length: 64" l (160 cm)
 Tabletop Material: Carbon fiber with integral head section
 Mattress: 2" thick (5 cm)
 Attenuation:
 Tabletop: .7 mm Aluminum equivalence
 Mattress: .5 mm Aluminum equivalence

Motions:

Height Adjustable: 30" to 40" (75 to 100 cm)
 Trendelenburg: 0° to 20°
 Reverse Trendelenburg: 0° to 20°
 Lateral Roll: 0° to 20°

Controls:

Hand Control: activates height, lateral roll and Trendelenburg motions

Optional:

Foot Control: activates height, lateral roll and Trendelenburg motions

Wheels:

Head End – 3" (7.62 cm) swivel casters, integral locking system on base
 Foot End – 5" (12.7 cm) swivel casters with central locking
 Accessory Rails: Standard OR accessory rails 30" x 1.12" x .375"
 (75 x 2.86 x .95 cm) mounted near foot end of table
 Patient Restraints: Three body straps
 Finish: Powder coat
 Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
 Lift Capacity: May vary between 460 lb - 500 lb (209 - 227 kg) based on line voltage
 Power: 115 VAC or 230 VAC
 Warranty: one year
 Certifications: ETL and cETL Listed to UL-2601-1®, CAN/CSA® C22.2
 No.: 601-1-M90, EN 60601-1 and IEC 60601-1-2

058-870	Table, Pain Management C-Arm - 870, 115 VAC	\$16,900.00
058-875	Table, Pain Management C-Arm - 870, 230 VAC	16,900.00

Optional:

058-877	Controller, Foot	\$895.00
056-850	IV Pole, Rail Mounted	390.00
056-865	Arm Board, Carbon Fiber, one arm	600.00
056-866	Arm Board, Carbon Fiber, two arms.....	1,150.00

Surgical C-Arm Table - 846

Affordably priced float-top table



SPECIFICATIONS:

Dimensions: 104" l x 26" w (260 x 66 cm) with OR accessory rails
 Tabletop: 101" l x 24" w (252 x 61 cm)
 Radiolucent Area:
 51" l x 24" w (128 x 61 cm) unobstructed head to accessory rail
 20" l x 17.5" w (50.8 x 44 cm) additional, between accessory rails
 Radiolucent Area, Total Length: 71" l (178 cm)
 Tabletop Material: Carbon Fiber with integral head section
 Mattress: 2" thick (5 cm)
 Attenuation:
 Tabletop: .7 mm Aluminum equivalence
 Mattress: .5 mm Aluminum equivalence

Motions:

Height: Adjustable from 33.5" to 43.5" (84 to 109 cm)
 Tabletop X Motion (head-to-toe float): 20" (51 cm)
 Tabletop Y Motion (side-to-side float): 10" (25.4 cm)

Controls:

Hand Control: height
 Push Button: releases X-Y Tabletop motions

Wheels:

Head End – 3" (7.62 cm) swivel casters, integral locking system on base
 Foot End – 5" (12.7 cm) swivel casters with central locking
 Accessory Rails: Standard OR accessory rails 49.5" x 1.12" x .375"
 (124 x 2.86 x .95 cm) mounted near foot end of table
 Patient Restraints: Three body straps
 Patient Capacity: 500 lb (227 kg); weight tested to four times the patient load rating.
 Lift Capacity: May vary between 460 lb - 500 lb (209 - 227 kg) based on line voltage
 Power: 115 VAC or 230 VAC
 Warranty: one year
 Certifications: ETL and cETL Listed to UL-2601-1®, CAN/CSA® C22.2
 No.: 601-1-M90 and EN 60601-1

058-846	Table, Surgical C-Arm - 846, 115 VAC	\$18,900.00
058-847	Table, Surgical C-Arm - 846, 230 VAC	18,900.00

Optional:

058-848	Catheter Tray Extension (removable).....	\$400.00
056-850	IV Pole, Rail Mounted	390.00
056-865	Arm Board, Carbon Fiber, one arm	600.00
056-866	Arm Board, Carbon Fiber, two arms.....	1,150.00